

Esco MIRI® Laminar Flow Cabinet

A vertical laminar flow cabinet that provides a constant and reliable airflow giving superior quality on the air supplying the working zone. Thereupon provides sample protection thoroughly. The MIRI® Laminar Flow is equipped with a highly energy-efficient ECM blower which minimize power consumption without affecting the cabinet's performance.



Key Features



ISO Class 3

Equipped with HEPA/ULPA filter with efficiency of 99.999% towards micro particulate of 0.3 µm.



Energy Efficient

The latest generation ECM blower allows the cabinet to have minimal power consumption while maintaining the cabinet's performance.



Anti-Microbial Coating

The ISOCIDE™ effectively inhibits and eliminates the growth of unwanted microbes at the surface of the cabinet's exterior.



Microscope Provision

The microscope provision allows user to easily install a stereo-microscope unit at their convenience.



Built-in Monitor

User may opt to have a touchscreen monitor equipped with a mini PC for better workflow.

Optional Accessories



Support Stand:

- With leveling feet
- With Caster Wheels
- Motorized stand with Caster Wheels



UV Kit



Want to know more about this product? Scan the QR Code

Esco Medical Aps
Kringelled 10, 8250 Egaa, Denmark
medical@escolifesciences.com
www.esco-medical.com

Esco Micro Pte Ltd (Headquarters)
19 Changi South Street 1, Singapore
486 779 • Tel +65 6542 0833

Esco Global Offices: Bangladesh | China | Denmark | Germany | Hong Kong | Indonesia | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam